

This Page Is Inserted by IFW Operations  
and is not a part of the Official Record

## **BEST AVAILABLE IMAGES**

Defective images within this document are accurate representations of the original documents submitted by the applicant.

Defects in the images may include (but are not limited to):

- BLACK BORDERS
- TEXT CUT OFF AT TOP, BOTTOM OR SIDES
- FADED TEXT
- ILLEGIBLE TEXT
- SKEWED/SLANTED IMAGES
- COLORED PHOTOS
- BLACK OR VERY BLACK AND WHITE DARK PHOTOS
- GRAY SCALE DOCUMENTS

**IMAGES ARE BEST AVAILABLE COPY.**

**As rescanning documents *will not* correct images,  
please do not report the images to the  
Image Problem Mailbox.**



Docket No. CDS-226

**IN THE UNITED STATES PATENT AND TRADEMARK OFFICE**

Applicants : Thomas J. Cummins, et al.

Serial No. : 09/675,828

Art Unit: 1637

Filed : September 29, 2000

Examiner: T.E.Strzelecka

For : DIAGNOSTIC COMPOSITIONS, ELEMENTS, METHODS AND  
TEST KITS FOR AMPLIFICATION AND DETECTION OF  
TWO OR MORE DNA's USING PRIMERS HAVING MATCHED  
MELTING TEMPERATURES

I hereby certify that this correspondence is being deposited with the  
United States Postal Service as first class mail in an envelope addressed  
to: Commissioner for Patents, P.O. Box 1450, Alexandria, VA 22313-1450 on

June 29, 2004

(Date)

Catherine Kurtz Gowen

Name of applicant, assignee, or Registered Representative

C. K. Gowen  
(Signature)

June 24, 2004  
(Date of Signature)

Commissioner for Patents  
P.O. Box 1450  
Alexandria, VA 22313-1450

**DECLARATION OF SUSAN MELISSA WERNER**  
**PURSUANT TO 37 CFR 1.182**

Dear Sir:

SUSAN M. WERNER DECLARES AS FOLLOWS:

1. I, SUSAN M. WERNER, am Regulatory Affairs Manager at  
Ortho-Clinical Diagnostics, Inc. and a co-inventor of the  
claims of the captioned application.

2. On the date the original Declaration was signed by me on May 14, 1993, I was known as Susan M. Atwood.

3. I was married on September 2, 2000, and have thereupon changed my name to Susan M. Werner.

4. Attorney of record Catherine Kurtz Gowen was first notified of my name change on or about May 25, 2004.


I hereby declare that all statements made herein of my own knowledge are true and that all statements made on information and belief are believed to be true; and further that these statements were made with the knowledge that willful false statements and the like so made are punishable by fine or imprisonment, or both, under Section 1001 of Title 18 of the United States Code and that such willful false statements may jeopardize the validity of the application or any patent issued thereon.

Susan M. Atwood  
Susan M. Atwood

Susan M. Werner  
Susan M. Werner

6/24/04  
Date

6/24/04  
Date



Docket No. CDS-226

IN THE UNITED STATES PATENT AND TRADEMARK OFFICE

Applicants : Thomas J. Cummins, et al.

Serial No. : 09/675,828

Art Unit: 1637

Filed : September 29, 2000

Examiner: T.E.Strzelecka

For : DIAGNOSTIC COMPOSITIONS, ELEMENTS, METHODS AND  
TEST KITS FOR AMPLIFICATION AND DETECTION OF  
TWO OR MORE DNA's USING PRIMERS HAVING MATCHED  
MELTING TEMPERATURES

I hereby certify that this correspondence is being deposited with the  
United States Postal Service as first class mail in an envelope addressed  
to: Commissioner for Patents, P.O. Box 1450, Alexandria, VA 22313-1450 on

June 29, 2004  
(Date)

Catherine Kurtz Gowen

Name of applicant, assignee, or Registered Representative

C. Gowen  
(Signature)

June 29, 2004  
(Date of Signature)

Commissioner for Patents  
P.O. Box 1450  
Alexandria, VA 22313-1450

DECLARATION UNDER 37 CFR 1.131

Dear Sir:

I, Susan M. Werner, depose and state:

1. I am coinventor of claims 28, 29, 33 and 36-40 of the above-identified patent application.

2. Prior to December 20, 1990, in collaboration with co-inventors Thomas J. Cummins, Lynn Bergemeyer, John

Bruce Findlay, John W.H. Sutherland and JoAnne H. Kerschner, I conceived and reduced to practice in this country, methods for the simultaneous amplification and detection for 2 or more target DNAs as described and claimed in our application. The following EXHIBITS A, B, C and D evidence the claimed invention:

A. A copy of my Laboratory Notebook No. BB1001, pages 89-90 and pages 95-97, (dated before December 20, 1990), provided as Exhibit A. The notebook entry shows the amplification and detection of 2 targets separated from each other along opposing strands by from 90-400 nucleotides, wherein the amplification uses 4 primers having Tms between 65°-74°C or 67°-74°C, wherein the primers are within about 2°C of each other, wherein priming and primer extension take place at 70°C (all the above documented specifically at pages 89-90 and 95-97) and wherein the primers have nucleotide lengths which differ from each other by no more than 5 nucleotides.

B. A copy of my Laboratory Notebook No. BB1001, pages 179-180 (dated after December 20, 1990 but before July 1992), provided as Exhibit B. Confirmation that the primers nucleotide lengths differed from each other by no more than 5 nucleotides is found at pages 179-180.

C. A copy of Kodak Technical Report (hereafter "TR") RPT. ACC. NO. 249038H co-authored by coinventor John W.H. Sutherland, (dated before December 20, 1990), title page and page 18, provided as Exhibit C. The title page provides an Abstract which summarizes the invention with respect to primer Tms and provides a working formula for calculating primer Tm as:

$$T_m(^{\circ}C) = 66.5 + 0.36 (\%G+C) - 384/N,$$

which calculation was used for obtaining primer Tms for the work done as documented in EXHIBIT A. EXHIBIT C goes on to state at page 18 that following the additional research described in the text of this TR, a modified regression equation for calculating primer Tm was derived as:

$$T_{m \text{ calc}} = 67.5 + 0.34 (\%G+C) - 395/N,$$

which it is said at the bottom of page 18 "differed only slightly from" the prior equation. This latter formula is identical to the formula provided in the patent disclosure at page 12, line 30 and page 26, line 30, as:

$$T_m(^{\circ}C) = 67.5 + 0.34(G+C) - 395/N,$$

and which formula appears in the pending claims.

D. A copy of co-inventor Thomas J. Cummins's Laboratory Notebook No. AA9223 at pages 13-14 (dated before December 20, 1990) provided as EXHIBIT D. The entry on these pages is a record of the conception of use of matched primer Tms in use of PCR that employs combined annealing and extension temperatures.


3. The date deleted from Exhibit A is prior to December 20, 2000 and July 1992.

4. The date deleted from Exhibit B is prior to July 1992.

5. The date deleted from Exhibit C is prior to December 20, 2000 and July 1992.

6. The date deleted from Exhibit D is prior to December 20, 2000 and July 1992.

I hereby declare that all statements made herein of my own knowledge are true and that all statements made on information and belief are believed to be true; and further that these statements were made with the knowledge that willful false statements and the like so made are punishable by fine or imprisonment, or both, under Section 1001 of Title 18 of the United States code and that such willful false statements may jeopardize the validity of the application or any patent issued thereon.

  
Susan M. Werner

Date 6/24/04

## RESEARCH / DEVELOPMENT

EASTMAN KODAK COMPANY

Date.

Problem:

HIV PCR

Objective: To determine if new batch of dNTP's functions under HIV assay conditions.

<b>REAGENTS:</b>			
10X HI SALT PCR BUFFER - pH 8.0 (100mM Mg++ , 100mM Tris) batch 006 prepared			
dNTP'S - Sigma batch 006 100µM prepared		by NG	
PRIMER 1 - EKGAG-F AA9316/106 10µM			
PRIMER 2 - EKGAG-R AA9316/106 10µM			
PRIMER 3 - EKENV-F AA9316/107 10µM			
PRIMER 4 - EKENV-R AA9316/120 10µM			
KODAK TAQ - (4 units/µl) EK004 prepared		by SA	
GENOMIC - H9 and 8E5 cell lysate #0018E5			
REAGENT	AMT/TUBE (µl)	# OF TUBES	STOCK AMT (µl)
10X PCR BUFFER	25	9	225
dNTP's	15	9	135
PRIMER 1	25	9	225
PRIMER 2	25	9	225
PRIMER 3	25	9	225
PRIMER 4	25	9	225
GLYCEROL (7.5%)	18.75	9	168.75
KTAQ (4 units/µl)	10	9	90
GENOMIC	25	9	225
TE BUFFER	56.25	9	506.25
TOTAL			2250
FILL POUCHES WITH 250µl			

Processor #4 profile: 90°C preheat - 10 seconds  
 96°C for 30 seconds  
 70°C for 60 seconds  
 for 40 cycles

~~Blank pouches #194~~

Signature

Susan M. Atwood

The foregoing disclosed to me on

Witness

RESEARCH / DEVELOPMENT

Notebook No.

EASTMAN KODAK COMPANY

BB1001

Date \_\_\_\_\_

Problem: page 89 cont'd


8 blank pouches were run:

<u>position</u>	<u>serial #</u>	<u>position</u>	<u>serial #</u>
1	194974	5	194964
2	194966	6	194973
3	194965	7	194968
4	194975	8	194971

Results

sample #

1 2 3 4 5 6 7 8



Both FK6AC and FKENV products amplified well;  
dNTP's don't seem to be a problem in the  
HIV system.

Signature

Susan M. Atwood

The foregoing disclosed to me on

[Signature]  
Witness



BB1001

EASTMAN KODAK COMPANY

Date.

Problem:

Day 0 Keeping

Began keeping study of tag polymerase in PCR mix solution. Using @ glycerol and HIV primers. Keeping at 3 different temperatures will be monitored; -20°C, 4°C and room temp. Tag activity assay will be performed by Dave Sharkey on the same day that I test for PCR ability.

10X HI SALT PCR BUFFER - pH 8.0 (100mM Mg++, 100mM Tris) batch 006 prepared

Stock

dNTP'S - Sigma 100mM prepared

by LB

Mix

PRIMER 1 - EKGAG-F oligo #699 10µM

made:

PRIMER 2 - EKGAG-R oligo #776 10µM

PRIMER 3 - EKENY-F oligo #700 10µM

PRIMER 4 - EKENY-R oligo #775 10µM

KODAK TAQ - batch 3B rebuild prepare

by SA

-20°C =

freezer

at U of R

4°C =

fridge

at U of R

Room temp =

desk in

room 5-8128

at U of R

REAGENT	AMT/TUBE (µl)	# OF TUBES	STOCK AMT (µl)
10X PCR BUFFER	25	52	1300
dNTP's	15	52	780
PRIMER 1	25	52	1300
PRIMER 2	25	52	1300
PRIMER 3	25	52	1300
PRIMER 4	25	52	1300
GLYCEROL (7.5%)	18.75	52	975
TE BUFFER	56.25	52	2925
TOTAL			11180

ALIUQUOT 215µL INTO THREE TUBES. TO EACH OF THESE TUBES ADD 10µL OF 40 UNITS/µL TAQ GIVEN TO DAVE SHARKEY AND STORED ONE TUBE EACH AT -20°C, 4°C, AND ROOM TEMP.

TO THE REMAINING MIX ADD 490µL OF 4 UNITS/µL TAQ. ALIUQUOT 225µL INTO 48 TUBES.

3 TUBES GIVEN TO DAVE SHARKEY. STORE 8 TUBES AT ROOM TEMP, 17 TUBES AT 4°C, AND 18 TUBES AT -20°C. BEFORE RUNNING PCR, ADD 25µL OF TARGET TO EACH TUBE.

TEST DATES:

DAY 0 = 11/9/90

1 WEEK = 11/16/90

1 MONTH = 12/10/90

2 MONTHS = 1/9/91

3 MONTHS = 2/8/91

4 MONTHS = 3/15/91

6 MONTHS = 4/15/91

1ul of 2.26% thimerisol was added to each tube to be stored at room temp - prevents bacterial growth.

Susan M. Atwood

*[Signature]*  
Witness

## RESEARCH / DEVELOPMENT

Notebook No.

EASTMAN KODAK COMPANY

BB1001

Date: \_\_\_\_\_

Problem: \_\_\_\_\_

Day 0 Keeping

Initial tag activity assay:

smp	raw cpm	cpm-bkg	dpm	nmol dA	U/ml	U/ml*dil'n	smp	U/ml	Tot U
1	302.5	41.50	189	0.000	0.043	0.04	No Enzyme	0.00	
2	8240.0	7979.00	36268	0.082	8.258	165.17			
3	8172.5	7911.50	35961	0.081	8.188	163.77	22°C/40U	162.85	36.8
4	8795.0	8534.00	38791	0.087	8.833	176.66	4°C/40U	169.04	38.0
5	10466.5	10205.50	46389	0.104	10.563	2112.56	-20°C/40U	176.50	39.7
6	11428.5	11167.50	50761	0.114	11.558	2311.70			
7	11320.0	11059.00	50268	0.113	11.446	2289.24	22°C/400U	2211.56	499.8
8	219.5	-41.50	-189	0.000	-0.043	-0.04	4°C/400U	2373.96	534.1
9	8016.0	7755.00	35250	0.079	8.027	160.53	-20°C/400U	2296.33	516.7
10	8682.0	8421.00	38277	0.086	8.716	174.32			
11	8780.0	8519.00	38723	0.087	8.817	176.35			
12	11423.0	11162.00	50736	0.114	11.553	2310.56			
13	12030.0	11769.00	53495	0.120	12.181	2436.21			
14	11388.5	11127.50	50580	0.114	11.517	2303.42			
15	23.0						CPI1	DPI1	uCi
16	28.0	bkg=	261.000			Liq =	16547328.0	42429046.2	13.82
17	977.5	=	samples 1,8			Flt =	118907.0	540486.4	0.18
18	4481.5					Tot =	16666235.0	42969532.5	
							238089.1		
							53.6		

Raw this initial mix made right away in 4 blast pouches. Added 25 ul of #0018E5 target to each tube, then loaded them into blast pouches.

Pouches were run on Processor #4 with the following profile: preheat = 90°C for 10 sec.  
then 40 cycles of  
96°C for 30 sec.  
70°C for 60 sec

target placed at -20°C for storage

Signature

Susan M. Atwood

EASTMAN KODAK COMPANY

Date

BB1001

Problem:

Day 0 Keeping

Pouches were loaded as follows:

position

Serial #

2

195074

(1)

3

195073\*

(2)

\* = PCR burst seal failure

6

195075

(3)

leakage during processing

7

195076


(4)

Results

Contents of PCR blister were removed into a tube and run

on an EtBr stained 4% agarose gel

3ul tracking dye + 12 ul sample → load well w/ 12 ul

Marker =  $\phi$ X174 Hae III  
digestSurecell Results

sample	AA7565-141 GAG	AA7565-142 ENV	AA7565-143 LTR	beads at 0.25%
1	6.5/6	8/7.5	0/0	
2	7/6	8.5/8	0/0	
3	7/8.5	7.5/9	0/0	
4	7.5/7.5	8.5/8.5	0/0	

Signature

Susan M. Howard

The foregoing disclosed to me on:

Witness

Notebook No.

## RESEARCH / DEVELOPMENT

EASTMAN KODAK COMPANY

Date: \_\_\_\_\_

BB1001

Problem: \_\_\_\_\_

HIV Oligo Preparation

Obtained primers from L. Kunzelbauer - diluted to 10 $\mu$ M  
in batch 006 TE buffer Aliquotted into 400  $\mu$ l volumes  
and stored at -20°C

CODE NAME: EKGAG-F #976 TRITYL: OFF  
PURPOSE: HIV primer  
SEQUENCE: AGT GGG GGG ACA TCA AGC AAGC CAT GCA A  
U= ~~5'ADT~~

LENGTH: 28

BASE COMP:

MW:(grams)

TOTAL MW:

A= 9

331.2

2980.8

C= 6

307.2

1843.2

G= 10

347.2

3472

T= 3

322.2

966.6

9262.6

O.D: 20 (1 OD=33ug)

ug DNA: 660

pmoles DNA: 71,254.29

 $T_m = 73.36$ 

CODE NAME: EKGAG-R

TRITYL:

Biotin

[STOCK] $\mu$ M: 10

PURPOSE: HIV primer

SEQUENCE: TTC CTG CTA TGT CAC TTC CCC TTG GTT C

U= 5'ADT

RESUSPEND IN: 7,125.43  
( $\mu$ l)

LENGTH: 28

BASE COMP:

MW:(grams)

TOTAL MW:

A= 2

331.2

662.4

C= 10

307.2

3072

G= 4

347.2

1388.8

T= 12

322.2

3866.4

8989.6

O.D: 25 (1 OD=33ug)

ug DNA: 825

pmoles DNA: 91,772.72

[STOCK] $\mu$ M: 10 $T_m = 70.79$ RESUSPEND IN: 9,177.27  
( $\mu$ l)

The foregoing disclosed to me on.

Supra M. Atwood  
Witness

## RESEARCH / DEVELOPMENT

Notebook No.

Date \_\_\_\_\_

EASTMAN KODAK COMPANY

BB 1001

Problem: \_\_\_\_\_

HIV Oligo Preparation

CODE NAME: EKENV-F # 977

TRITYL: OFF

PURPOSE: HIV primer

SEQUENCE: GGT TCT TGG GAG CAG CAG GAA GCA CTA T

U= ~~5'ADTS~~

LENGTH: 28

BASE COMP:

MW:(grams)

TOTAL MW:

A= 7

331.2

2318.4

C= 5

307.2

1536

G= 10

347.2

3472

T= 6

322.2

1933.2

9259.6

O.D: 20 (1 OD=33ug)

ug DNA: 660

pmoles DNA: 71,277.38

 $T_m = 72.07$ 

[STOCK]uM: 10

RESUSPEND IN: 7,127.74  
(ul)

CODE NAME: EKENV-R # 948

TRITYL: BIOTIN

PURPOSE: HIV primer

SEQUENCE: CTT GAT GCC CCA GAC TGT GAG TTG CAA C

U= 5'ADTS

LENGTH: 28

BASE COMP:

MW:(grams)

TOTAL MW:

A= 6

331.2

1987.2

C= 8

307.2

2457.6

G= 7

347.2

2430.4

T= 7

322.2

2255.4

9130.6

O.D: 20 (1 OD=33ug)

ug DNA: 660

pmoles DNA: 72,284.41

[STOCK]uM: 10

 $T_m = 72.07$ RESUSPEND IN: 7,228.44  
(ul)

Signature \_\_\_\_\_

Susan M. Wood  
Debra J. Jure  
Witness



CONFIDENTIAL INFORMATION  
CONTROLLED DISTRIBUTION  
RECLASSIFY \_\_\_\_\_

RPT. ACC. NO.  
249038H

DATE

THIS COPY FOR:

# TECHNICAL REPORT

DIVISION  
Life Sciences

KODAK RESEARCH LABORATORIES  
ROCHESTER, N.Y.

STUDIES OF DNA DENATURATION AND RENATURATION IN POLYMERASE CHAIN REACTION BUFFER: PART V.  $T_m$ 'S OF COMMONLY USED PRIMERS, AND RELATIONSHIP TO AMPLIFICATION YIELDS IN PCR

AUTHORS

A. R. Mack and J.W.H. Sutherland

CONTRIBUTORS

Dr. John Findlay and  
Mrs. Marlene King

## ABSTRACT

$T_m$ 's in PCR buffer have been measured for a variety of primers and probes currently used on the AIDS project. The  $T_m$ 's for the AIDS primers SK38 and 39 (66°C and 65.5°C) are not very different from those for beta-globin primers PC03 and PC04 ( $T_m$ 's = 64.5 and 67°C, respectively).

Using a basis set of these 11 oligomers for  $T_m$  calculations, the following theoretical relationship between % guanine and cytosine content as well as number of base pairs,  $N$ , was determined:

$$T_m(^{\circ}\text{C}) = 66.5 + 0.36 (\%G + C) - 384/N$$

Although similar in functional form, this relationship differs quantitatively from analogous expressions reported in the literature.

Actual melting transitions generally were found to be in reasonable agreement with an all-or-none (i.e., thermodynamic) model using two adjustable parameters, viz., the enthalpy and entropy changes of melting. This suggests that at any given temperature these short oligonucleotides and their complements existed in equilibrium between melted and annealed forms. Hence, partially melted species did not appear to exist in detectable amounts.

Most of the primers examined here have  $T_m$ 's which are less than the extension temperature of 70°C. When the yields obtained using these primers in PCR were supplemented with the  $T_m$ 's measured here, it became apparent that good PCR yields required annealing at temperatures at or below primer  $T_m$ 's. Respectable yields could be achieved with extension temperatures at or above  $T_m$ , but in this case most of the extension probably took place during the annealing stage. This suggests that we should either (1) lengthen the primers to raise their  $T_m$ 's above the extension temperature or (2) perform two-temperature PCR's, with the lower temperature at or below the  $T_m$  of the existing (unlengthened) primers.

APPENDIX - A Test of The Regression Equation Using  
Additional Oligonucleotides

RPT. ACC. NO. 249038H

CONTROLLED DISTRIBUTION

## APPENDIX - A Test of The Regression Equation Using Additional Oligonucleotides

Following the completion of the work described in the previous text, the regression equation obtained therein (Eq. 1) was used to predict the  $T_m$  of 10 additional oligonucleotides which were not contained in the original basis set described in Table 1. The length, composition, and measured  $T_m$ 's of the additional oligonucleotides is included in Table A2. The residuals ( $T_m \text{ calc} - T_m \text{ exptl}$ ) are given in Table A1(below):

Table A1: A Test of The Regression Equation (1)

oligonucleotide	Residual ( $^{\circ}\text{C}$ )
ARM-8	0.7
ARM-9	1.0
ARM-10	0.0
ARM-11	-1.0
ARM-12	1.1
ARM-6	-0.4
IBI	-2.1
FOBS-3	-0.5
PC03-4	1.8
PC04-4	-2.5

It can be seen that the predictions are generally fairly good (the average residual being only  $-0.2^{\circ}\text{C}$ ).

This data was then pooled with that of Table 1, generating a pooled set of 21 oligonucleotides. The modified regression equation for this new set:

$$T_m \text{ calc} = 67.5 + 0.34(\%G+C) - 395/N \quad (\text{A-1})$$

differed only slightly from Eq. (1):

$$T_m \text{ calc} = 66.5 + 0.36(\%G+C) - 384/N$$

Residuals and fits are shown in Table A-2 and Figs. A-1 and A-2.

Notebook No.

## RESEARCH / DEVELOPMENT

EASTMAN KODAK COMPANY

Date \_\_\_\_\_

Problem:

*Coamplification of Multiple Primers*

Send to: Dan  
From: Thomas Cummins  
Subject: Patents

These are important patent items:

1. Common Tms for coamplification of two or more primer sets using PCR.

Outline. When applying two temperature PCR whereby the annealing and extension temperature are combined, it is important that multiple primer sets are structured so that the Tm of the primers are equivalent. Tm is dictated primarily by primer length and G / C content. Equal Tms (w/in 1 - 3 degrees) will ensure that multiple primer pairs anneal at roughly the same efficiency so that the overall efficiency of the amplification reaction is not altered for any one of the multiple primers. This ensures that the sensitivity of the primers is equal.

advantages....Fast PCR = better specificity, faster test time. Common Tms for coamplified primers = allows multiple sets of primers to be coamplified in a single reaction using two temperature PCR ensuring equivalent amplification efficiencies for all coamplified primers.

efficiencies for all coamplified primers.

Example:

Targets HIV gag gene and env gene:

Primers = gag primer = Tm 68 C..... env primer = Tm = 58 C...if these primer sets are coamplified using 2 temp PCR, where the annealing / extension temperature is 67 C, only the gag primer will amplify because the temperature is well above the melting point of the env primer and the primer will not bind at all to the target sequence or bind at a much lower efficiency than the gag primer leading to a compromised amplification.

By bringin the env Tm to 68 C (similar to the gag primer) both sets can now be amplified at the same efficiency at the single temperature of 67 C.

Thanks, Tom

18226-5/86 L.P.S.

Signature

The foregoing disclosed to me on.

Witness



RESEARCH / DEVELOPMENT  
EASTMAN KODAK COMPANY

Notebook No.

Date

AA9223

AA9223

Problem:

*Similar Tms for the detection of amplified*

SEND A NOTE

E04

Send to: Dan  
From: Thomas Cummins  
Subject: Patent

2. Common probe Tms for the codetection of two or more amplified nucleic acid sequences following the polymerase chain reaction.

Outline: When detection two or more amplified nucleic acid sequences (products) following coamplification using the polymerase chain reaction, it is an advantage to have the probe/product complex for each primer set constructed so that the Tms are roughly the same (w/in 1 - 3 degrees). Similar Tms will thus allow a single wash stringency (temperature and ionic strength) to be used during the detection chemistry. This is desirable in a system where a single chemistry is to be used for the detection of multiple products.

Thanks, Tom

Signature

The foregoing disclosed to me on.

Witness